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
2. ASTM Standard F 136, 2013, "Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)," ASTM International, West Conshohocken, PA, 2013, DOI: 10.1520/F0136, www.astm.org.

3. Zamanian, K., & van den Dolder, I. (2014, May 17). Clinical Implantology- Conical internal connections will fuel growth in dental implant market. Retrieved August 1, 2014, from Dental Tribune: www.dental-tribune.com

4. BioHorizons' Internal Reference for Manufactured Components per Year and Fractured Implants with Internal Conical Connection since 2010.

5. Wright, R. B. (2005, June). How Tight Should Fasteners be Tightened? Fastner Technology International , pp. 24-25.

6. Underwood, J. R., & Chiuiini, M. (1998). Bolted Connection. In Structural Design: A Practical Guide for Architects (pp. 223-225). New York: John Wiley & Sons.

7. Aquirrebeitia, J., & Abasolo, M. (2013). Dental Implants with Conical Implant-Abutment Interface: Influence of the Conical Angled Difference on the Mechanical Behavior of the Implant. Int J Oral Maxillofac Implants.

8. Bozkaya, D., & Muftu, S. (2005). Mechanics of the taper integrated screwed-in (TIS) abutments used in dental implants. Journal of Biomechanics , 38:87-97.

9.Alexander, H., Ricci, J. L., & Hrico, G. J. (2009). Mechanical Basis for Bone Retention around Dental Implants. J Biomed Mater Res Part B: Appl Biomater , 306-311.

10.Sarmah A, Mathur AK, Gupta V, Pai VS, Nandini S. Finite Element Analysis of Dental Implant as Orthodontic Anchorage. J Contemp Dent Pract 2011; 12(4):259-264

11.Soncini, M. (2003). Experimental Tests for Evaluating Ultimate Load. In A. N. Natali, Dental Biomechanics (p. 116). London: Taylor & Francis.

12. Seabra Barbosa, G. A., Bernardes, S. R., Domingues das Neves, F., Fernandez Neto, A. J., Chiarello de Mattos, M. d., & Ribeiro, R. F. (2008). Relation between Implant/ Abutment Vertical Misfit and Torque Loss of Abutment Screws. Braz. Dent J. , vol 19 No 4.

13. Coppede, A. R., Chiarello de Mattos, M. d., & Aparecida de Mattias Sartori, I. (2009). Fracture resistance of the implant-abutment connection in implants with internal hex and internal conical connections under oblique compressive loading: an in vitro study. International Journal of Prosthodontics , 22:283-286.

14. Ugurel, C. S., Steiner M., Isik-Ozkol G., Kutay O., & Kern M. (2013). Mechanical Resistance of Screwless Morse Taper and Screw-Retained Implant-Abutment Connections. Clin Oral Imp Research. 10.1111 /clr.12303

15. Shigley, J. E. (1989). The design of screws, fastners and connections. In Mechanical Engineering Design (p. 337). McGraw-Hill Inc.

16. Huei-Ming Chai (2005). Application of Stress-Strain Curve. Measurement of Stress and Strain. www.pt.ntu.edu.tw.

17. Spiralock. (2014). Retrieved July 24, 2014, from Stanley Engineered Fastening: http://www.stanleyengineeredfastening.com

18. Fatigue Testing for Endosseous Dental Implants. American National Standard/American Dental Association. Standard No. 127.Identical adoption of ISO 14801, 2007.

19. Class II Special Controls Guidance Documents: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments. U.S. Food and Drug Administration. Guidance for Industry and FDA Staff. May 2002.

20. Boggan, R. S., Strong, J. T., Misch, C. E., & Bidez, M. W. (1999). Influence of Hex Geometry and Prosthetic Table Width on Static and Fatigue Strength of Dental Implants. The Journal of Prosthetic Dentistry , 436-440.

21. Turzo, K. (2012). Surface Aspects of Titanium Dental Implants, Biotechnology - Molecular Studies and Novel Applications for Improved Quality of Human Life, Prof. Reda Sammour (Ed.), ISBN: 978-953-51-0151-2.

22. ASTM Standard F 67, 2013, "Standard Specification for Unalloyed Titanium, for Surgical Implant Applications (UNS R50250, UNS R50400, UNS R50550, UNS R50700)," ASTM International, West Conshohocken, PA, 2013, DOI: 10.1520/F0067, www.astm.org.

23. Hallab, N. J., Jacobs, J. J., & Katz, J. L. (2004). Orthopedic Applications. In Biomaterials Science: An introduction to Materials in Medicine (p. 537). California: Elsevier Academic Press.

24. Okazaki, Y. (2012). Comparison of Fatigue Properties and Fatigue Crack Growth Rates of Various Implantable Metals. Materials , 2982-3001.

25. Aran, A. (2007). Manufacturing Properties of Engineering Materials. ITU Department of Mechanical Engineering.

26. Dental Implants with an Internal Connection. (2007, Janurary). Retrieved August 4, 2014, from OsseoNews: www.osseonews.com

BioHorizons® and Laser-Lok® are registered trademarks of BioHorizons. Spiralock® is a registered trademark of Spiralock Corporation. Not all products shown or described in this literature are available in all countries. As applicable, BioHorizons products are cleared for sale in the European Union under the EU Medical Device Directive 93/42/EEC and the tissues and cells Directive 2004/23/EC. We are proud to be registered to ISO 13485:2003, the international quality management system standard for medical devices, which supports and maintains our product licences with Health Canada and in other markets around the globe. Original language is English. ©BioHorizons. All Rights Reserved.

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