

テーパードティッシュレベルインプラント

インプラント体直径	 3.0mm	 3.8mm	 4.6mm	 5.8mm	
補綴プラットフォーム直径	 3.5mm	 3.5mm	 4.5mm	 5.7mm	
Laser-Lokゾーン	2.0mm	2.0mm	2.0mm	2.0mm	
インプラント先端部径	2.0mm	2.8mm	3.1mm	3.9mm	
カラー部の高さ	2.3mm	2.3mm	2.3mm	2.3mm	
ベベル部の高さ	0.5mm	0.5mm	0.5mm	0.5mm	
カラー部の最大径	4.5mm	4.5mm	5.5mm	6.7mm	
インプラント体の長さ	7.5mm	—	TPLY3807	TTLG4607	TTLB5807
	9.0mm	—	TPLY3809	TTLG4609	TTLB5809
	10.5mm	TPLY3010	TPLY3810	TTLG4610	TTLB5810
	12.0mm	TPLY3012	TPLY3812	TTLG4612	TTLB5812

